

Before & After



Courtesy of A. Gaspar, MD



Courtesy of R. Valdivia Sing, MD

Get the eyebrows you've always wanted!

Impressive and natural-looking results

Three to four sessions with a three-week interval are recommended in order to achieve long-lasting effects.

110717

Fotona
choose perfection



VectorLift®

Non-invasive laser eyebrow lifting

What is VectorLift®?

Never underestimate the power of beautiful eyebrows to give your face a more youthful and attractive appearance. VectorLift® is an innovative solution for **non-invasive eyebrow lifting** and **upper eyelid rejuvenation**.

What is the treatment like?

VectorLift® is a comfortable and non-invasive procedure, allowing you to immediately return to your normal daily routine. The total treatment time is around 25 minutes, and it is usually performed under topical anesthesia. Some redness or light swelling may be present for at most 12 hours after the procedure, and some light skin peeling may occur.

Key advantages of VectorLift®

- Forehead eyebrow elevation
- Eyelid tightening
- Natural-looking results
- Non-invasive procedure
- No downtime



How does it work?

VectorLift® combines forehead eyebrow elevation with eyelid and scalp tightening. The first step involves treating the forehead, followed by extending the treatment to the scalp in order to provide a strong hold or anchor point for the lifting effect. Finally, treatment of the upper eyelid contributes to addressing eyelid laxity.



VectorLift® provides the optimal combination of an immediate skin tightening effect with long-term rejuvenation due to new collagen formation, all achieved using an innovative hyperstacking technique of low-energy Fotona SMOOTH® pulses, which are safe and gentle.